



# Recombinant K. pneumoniae Cardiolipin synthase (partial) (DAG-WT4167)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                  |                                                                                                                               |
|------------------|-------------------------------------------------------------------------------------------------------------------------------|
| Product Overview | This item requires custom production and lead time is between 4-12 weeks. We can custom produce based on your specifications. |
| Species          | K. pneumoniae                                                                                                                 |
| Purity           | > 85% as determined by SDS-PAGE                                                                                               |
| Conjugate        | TBD                                                                                                                           |
| Applications     | Immunoassays                                                                                                                  |
| Format           | Liquid or Lyophilized powder                                                                                                  |
| Size             | 1 mg                                                                                                                          |
| Buffer           | Tris/PBS-based buffer                                                                                                         |
| Storage          | Store at -20°C to -80°C                                                                                                       |

## BACKGROUND

**Keywords** Klebsiella pneumoniae; K. pneumoniae; Cardiolipin synthase