



Trifloxystrobin [BSA] (DAG-WT421)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Aluminum-ETDA conjugated with OVA
Specificity	This product is ELISA-paired with CABT-L7816
Species	N/A
Conjugate	BSA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Trifloxystrobin is the methyl ester of (2E)-(methoxyimino)[2-({[(E)-{1-[3-(trifluoromethyl)phenyl]ethylidene}amino]oxy}methyl)phenyl]acetic acid. A foliar applied fungicide for cereals which is particularly active against Ascomycetes, Deuteromycetes and Oomycetes It has a role as a mitochondrial cytochrome-bc1 complex inhibitor and an antifungal agrochemical. It is an oxime O-ether, an organofluorine compound, a methyl ester and a methoxyiminoacetate strobilurin antifungal agent.
Keywords	Trifloxystrobin; Trifloxystrobine; Fungicides; BSA conjugates