



Abietic acid [BSA] (DAG-WT427)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------|---|
| Specificity | This product is ELISA-paired with CABT-L7819 |
| Species | N/A |
| Conjugate | BSA |
| Format | Liquid |
| Concentration | Batch dependent - please inquire should you have specific requirements. |
| Size | 1 mg |
| Buffer | 0.01M PBS, pH7.4 |
| Preservative | None |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |

BACKGROUND

| | |
|--------------|---|
| Introduction | Abietic acid is an abietane diterpenoid that is abieta-7,13-diene substituted by a carboxy group at position 18. It has a role as a plant metabolite. It is an abietane diterpenoid and a monocarboxylic acid. It is a conjugate acid of an abietate. |
| Keywords | Abietic acid; Resin; Abietane diterpenoid; BSA conjugates |