



Cyclamate [BSA] (DAG-WT429)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This product is ELISA-paired with CABT-L7820
Species	N/A
Conjugate	BSA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Cyclohexylsulfamic acid is a member of the class of sulfamic acids that is sulfamic acid carrying an N-cyclohexyl substituent. It has a role as a human xenobiotic metabolite and an environmental contaminant. It derives from a sulfamic acid. It is a conjugate acid of a cyclohexylsulfamate.
Keywords	Cyclamate; Sweetener; Sodium cyclamate; BSA conjugates