



## Vitamin K3 [BSA] (DAG-WT431)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Specificity	This product is ELISA-paired with CABT-L7823
Species	N/A
Conjugate	BSA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

Introduction	Menadione is a member of the class of 1,4-naphthoquinones that is 1,4-naphthoquinone which is substituted at position 2 by a methyl group. It is used as a nutritional supplement and for the treatment of hypoprothrombinemia. It has a role as a nutraceutical, a human urinary metabolite, an angiogenesis inhibitor, an EC 3.4.22.69 (SARS coronavirus main proteinase) inhibitor and an antineoplastic agent. It is a member of 1,4-naphthoquinones and a vitamin K.
Keywords	Vitamin K3; Vitamin K; Menadione; BSA conjugates