



β-lactam receptor [BSA] (DAG-WT435)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This product is ELISA-paired with CABT-L7825
Species	N/A
Conjugate	BSA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	β-Lactam antibiotics inhibit bacteria by binding covalently to PBPs in the cytoplasmic membrane.
Keywords	β-lactam receptor; Beta-lactam; Antibiotics; BSA conjugates