



## Procymidone [OVA] (DAG-WT450)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Specificity	This product is ELISA-paired with CABT-L7805
Species	N/A
Conjugate	OVA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

Introduction	Procymidone is a pesticide. It is often used for killing unwanted ferns and nettles, and as a dicarboximide fungicide for killing fungi, for example as seed dressing, pre-harvest spray or post-harvest dip of lupins, grapes, stone fruit, strawberries.
Keywords	Procymidone; Pesticide; OVA conjugates