



## Chlordimeform [BSA] (DAG-WT455)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Specificity</b>	This product is ELISA-paired with CABT-L7808
<b>Species</b>	N/A
<b>Conjugate</b>	BSA
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M PBS, pH7.4
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

<b>Introduction</b>	Chlordimeform is an acaricide (pesticide) active mainly against motile forms of mites and ticks and against eggs and early instars of some Lepidoptera insects.
<b>Keywords</b>	Chlordimeform; Acaricide-insecticide; BSA conjugates