



## Cyproconazole [BSA] (DAG-WT459)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Specificity	This product is ELISA-paired with CABT-L7810
Species	N/A
Conjugate	BSA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

Introduction	Cyproconazole is an agricultural fungicide of the class of azoles, used on cereal crops, coffee, sugar beet, fruit trees and grapes, on sod farms and golf courses and on wood as a preservative.
Keywords	Cyproconazole; Fungicide; BSA conjugates