



Myclobutanil [BSA] (DAG-WT471)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This product is ELISA-paired with CABT-L7793
Species	N/A
Conjugate	BSA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	<p>C₁₅H₁₇ClN₄; CAS Number: 88671-89-0.</p> <p>Myclobutanil is a light yellow solid used as a fungicide. It is a triazole chemical used as a fungicide. It is a steroid demethylation inhibitor, specifically inhibiting ergosterol biosynthesis.</p>
Keywords	Myclobutanil; Fungicide; BSA conjugates