



## Triapenthenol [BSA] (DAG-WT479)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Specificity</b>	This product is ELISA-paired with CABT-L7797
<b>Species</b>	N/A
<b>Conjugate</b>	BSA
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M PBS, pH7.4
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

<b>Introduction</b>	Triapenthenol is a herbicide. It is highly soluble in water. Little is known regarding its environmental fate but it is considered to be persistent in soils. It is not highly toxic to mammals. It shows a high level of toxicity to birds. There is little information on its toxicity to biodiversity.
<b>Keywords</b>	Triapenthenol; Herbicide; BSA conjugates