



# Recombinant Human SHBG (a.a 30-402) [His] (DAG-WT514)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	Human
Purity	> 85% , as determined by SDS-PAGE
Conjugate	His
Applications	N/A
Molecular Weight	42 kDa
Format	Lyophilized Powder
Size	50 µg
Storage	Short Term: 2-8°C. Long Term: -20°C

## BACKGROUND

Introduction	Sex hormone-binding globulin (SHBG) or sex steroid-binding globulin (SSBG) is a glycoprotein that binds to androgens and estrogens. Other steroid hormones such as progesterone, cortisol, and other corticosteroids are bound by transcortin. SHBG is found in all vertebrates apart from birds.
Keywords	Sex hormone-binding globulin; SHBG; Steroid hormone

## GENE INFORMATION

UniProt ID	BC101785
------------	----------