



Recombinant Human IGF Binding Proteins-2 (IGFBP-2) (DAG-WT517)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% , as determined by HPLC
Conjugate	Unconjugated
Applications	N/A
Molecular Weight	32.4 kDa
Endotoxin	≤0.1 EU/μg
Format	Lyophilized Powder
Size	100 μg
Storage	Store at 2°C to 8°C.

BACKGROUND

Introduction	IGFBP-2 is a member of a highly conserved family of six insulin-like growth factor (IGF) binding proteins (IGFBPs) that modulate IGF-I actions, and they play an important role in the regulation of several cellular processes.
Keywords	Insulin-like growth factors; IGFBP2; Insulin