



# Recombinant Human IGFBP3 [His] (DAG-WT518)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

**Purity** > 95% as determined by SDS-PAGE

**Conjugate** His

**Applications** ELISA

**Molecular Weight** 29.6 kDa

**Format** Liquid

**Size** 50 µg, 500 µg, 1 mg

**Buffer** 40mM Sodium Citrate, pH 4.0

**Preservative** None

**Storage** Store at -20°C to -80°C

## BACKGROUND

**Keywords** IGFBP3; Insulin-like growth factor-binding protein 3

## GENE INFORMATION

**UniProt ID** P17936