



Recombinant Human IGF Binding Proteins-5 (IGFBP-5) (DAG-WT520)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% , as determined by SDS-PAGE
Conjugate	Unconjugated
Applications	N/A
Molecular Weight	28 kDa
Endotoxin	≤0.1 EU/μg
Format	Lyophilized Powder
Size	100 μg
Storage	Store at 2°C to 8°C.

BACKGROUND

Introduction	IGFBP-5 is a multifunctional protein that is capable of regulating IGF signaling both positively and negatively in different tissues and cells. It can also promote, or inhibit cell survival, proliferation, migration, etc. via mechanisms independent of IGF binding.
Keywords	Insulin-like growth factors; IGFBP5; Insulin