



# Recombinant Human TIM-3 [hFc] (DAG-WT5803)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	Human
Purity	> 90% as determined by SDS-PAGE
Conjugate	hFc
Applications	Immunoassays
Molecular Weight	46.6 kDa
Format	Lyophilized powder
Size	100 µg, 1 mg
Buffer	Lyophilized from sterile PBS, pH 7.4
Preservative	None
Storage	Store at -20°C to -80°C

## BACKGROUND

Keywords	TIM-3; HAVCR2; Hepatitis A virus cellular receptor 2
----------	--

## GENE INFORMATION

UniProt ID	NP_116171.3
------------	-------------