



Recombinant Lymphocytic choriomeningitis virus Glycoprotein G2 [His] (DAG-WT5841)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 85% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA
Molecular Weight	33.3 kDa
Format	Liquid or Lyophilized powder
Size	20 µg, 100 µg, 1 mg
Buffer	Tris/PBS-based buffer
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords Lymphocytic choriomeningitis virus; Glycoprotein G2; Pre-glycoprotein polyprotein GP complex

GENE INFORMATION

UniProt ID	P09991
------------	--------