



Recombinant Canine coronavirus Nucleoprotein (N) (DAG-WT6063)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 85% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	45.5 kDa
Format	Liquid or Lyophilized powder
Size	20 µg, 100 µg, 1 mg
Buffer	Tris/PBS-based buffer
Storage	Store at -20°C to -80°C. Avoid repeated freeze/thaw cycles

BACKGROUND

Keywords	Canine coronavirus; Nucleoprotein; N protein
-----------------	--

GENE INFORMATION

UniProt ID	Q7T6S8
-------------------	--------