



Recombinant SFTS Virus Nucleoprotein (DAG-WT640)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95% as determined by SDS-PAGE
Conjugate	TBD
Applications	ELISA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	0.01 M PBS, pH7.3
Preservative	None
Storage	Store at -20°C to -80°C

BACKGROUND

Introduction SFTSV is a Phlebovirus in the family of Bunyaviridae. It appears to be more closely related to the Uukuniemi virus serogroup than to the Sandfly fever group. It may be the prototype of a new serogroup within the Phlebovirus genus as it is only distantly related to both the existing serogroups.

Keywords SFTS Virus; SFTSV; Nucleoprotein