



MMAE [BSA] (DAG-WT669B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	10mM PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	MMAE (monomethyl auristatin E) is an antimitotic agent that binds to tubulin and inhibits the assembly of microtubules. Mc-VC-PAB-MMAE is a thiol reactive derivative of MMAE that containing a cathepsin cleavable valine-citruline p-aminobenzylcarbamate linkage. Mc-VC-PAB-MMAE is used to prepare antibody-drug conjugates, such as brentuximab vedotin.
--------------	---