



Vat-Cit-PAB-MMAD [KLH] (DAG-WT676K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Vat-Cit-PAB-MMAD conjugated with KLH
Conjugate	KLH
Cas.No	1415329-13-3
Applications	N/A
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	PBS buffer
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles

BACKGROUND

Introduction	MMAD is a tubulin inhibitor and used to produce site-specific antibody drug conjugates.
Keywords	Monomethyl auristatin D; MMAD; Auristatin