



Mc-MMAF [KLH] (DAG-WT678K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mc-MMAF conjugated with KLH
Conjugate	KLH
Cas.No	863971-19-1
Applications	N/A
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	PBS buffer
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles

BACKGROUND

Introduction	Monomethyl auristatin F (MMAF) is a synthetic antineoplastic agent. It is part of the approved drug belantamab mafodotin in multiple myeloma and some experimental anti-cancer antibody-drug conjugates such as vorsetuzumab mafodotin and SGN-CD19A. In International Nonproprietary Names for MMAF-antibody-conjugates, the name mafodotin refers to MMAF plus its attachment structure to the antibody.
Keywords	Monomethylauristatin F; Monomethyl auristatin F; MMAF