



FITC [BCP] (DAG-WT6832)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------|---------------------|
| Conjugate | BCP |
| Applications | Immunoassays |
| Format | Liquid |
| Size | 1 mg |
| Buffer | 10mM PBS, pH7.0±0.2 |
| Preservative | None |
| Storage | Store at 2°C to 8°C |

BACKGROUND

| | |
|-----------------|----------------------------------|
| Keywords | Fluorescein Isothiocyanate; FITC |
|-----------------|----------------------------------|