



DM1-SMCC [KLH] (DAG-WT683K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	DM1-SMCC conjugated with KLH
Conjugate	KLH
Cas.No	1228105-51-8
Applications	N/A
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	PBS buffer
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles

BACKGROUND

Introduction	DM1 (mertansine) is a thiol-containing derivative of maytansine. It binds to tubulin and inhibits the assembly of microtubules. DM1 is used to prepare antibody-drug conjugates, such as trastuzumab emtansine (abbreviated as T-DM1).
Keywords	Mertansine; DM1; Maytansine