



## PAC [BSA] (DAG-WT686B)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Product Overview</b>	PAC conjugated with BSA
<b>Conjugate</b>	BSA
<b>Cas.No</b>	2158322-33-7
<b>Applications</b>	N/A
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements
<b>Size</b>	1 mg
<b>Buffer</b>	PBS buffer
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles

### BACKGROUND

<b>Introduction</b>	PAC, consists the ADCs linker and PROTACs, conjugated to an antibody. PAC extracts from patent WO2017201449A1, compound LP2. PAC conjugated to an antibody is a more marked estrogen receptor-alpha (ERα) degrader compared to PROTAC (without Ab).
<b>Keywords</b>	PAC; PROTAC; ERα