



## PAC [KLH] (DAG-WT686K)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|                  |                                                                                      |
|------------------|--------------------------------------------------------------------------------------|
| Product Overview | PAC conjugated with KLH                                                              |
| Conjugate        | KLH                                                                                  |
| Cas.No           | 2158322-33-7                                                                         |
| Applications     | N/A                                                                                  |
| Format           | Liquid                                                                               |
| Concentration    | Batch dependent - please inquire should you have specific requirements               |
| Size             | 1 mg                                                                                 |
| Buffer           | PBS buffer                                                                           |
| Preservative     | None                                                                                 |
| Storage          | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles |

### BACKGROUND

|              |                                                                                                                                                                                                                                                              |
|--------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Introduction | PAC, consists the ADCs linker and PROTACs, conjugated to an antibody. PAC extracts from patent WO2017201449A1, compound LP2. PAC conjugated to an antibody is a more marked estrogen receptor-alpha (ER $\alpha$ ) degrader compared to PROTAC (without Ab). |
| Keywords     | PAC; PROTAC; ER $\alpha$                                                                                                                                                                                                                                     |