



Human Fibrinogen 3 (DAG-WT703)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Plasminogen, von Willebrand Factor and Fibronectin Depleted
Purity	> 95 % as determined by SDS-PAGE
Conjugate	N/A
Applications	N/A
Bio-activity	>95% clottable
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 g
Buffer	20mM sodium citrate-HCl buffer, pH 7.4
Preservative	None
Storage	Store at <-65°C. Avoid repeated freezing and thawing

BACKGROUND

Introduction	Fibrinogen (factor I) is a glycoprotein in vertebrates that helps in the formation of blood clots. It consists of a linear array of three nodules held together by a very thin thread which is estimated to have a diameter between 8 and 15 Angstrom (Å).
Keywords	Fibrinogen; Factor I