



Inactivated Natural West Nile Virus (Strain NY 2001-6263) Quality Control (DAG-WT7104)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	N/A
Applications	Molecular control, Immunoassays
Stability	24 months
Format	Liquid
Concentration	10 ⁵ –10 ⁷ TCID ₅₀ /mL 10 ⁷ –10 ¹⁰ copies/mL
Size	250 µL, 1 mL
Storage	Store at -80°C
Ship	Dry ice

BACKGROUND

Keywords	West nile virus; WNV
-----------------	----------------------