



Biotinylated Recombinant Human Neuropilin-1 [His,Avi] (DAG-WT713)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	A DNA sequence encoding human neuropilin-1 (a.a 22-644) was expressed with a polyhistidine tag and an Avi tag at the C-terminus.
Purity	> 95 % as determined by SDS-PAGE
Conjugate	His, Avi
Applications	ELISA
Molecular Weight	73.5 kDa
Format	Lyophilized
Size	25 µg
Buffer	Lyophilized from 0.22 µm filtered solution in PBS, pH7.4
Preservative	None
Storage	Store at -20°C to -80°C

BACKGROUND

Introduction	Neuropilin-1 is a protein that in humans is encoded by the NRP1 gene. In humans, the neuropilin 1 gene is located at 10p11.22. This is one of two human neuropilins.
Keywords	Neuropilin-1; NRP1

GENE INFORMATION

