



Inactivated Natural Salmonella spp. Quality Control (DAG-WT7198)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	N/A
Applications	Immunoassays
Stability	24 months
Format	Liquid
Concentration	10^6 – 10^9 CFU/mL
Size	250 μ L, 1 mL
Storage	Store at -80°C
Ship	Dry ice

BACKGROUND

Keywords	Salmonella; Salmonella spp.
-----------------	-----------------------------