



## Recombinant Human SLA/LP (DAG-WT7317)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Species</b>	Human
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Conjugate</b>	TBD
<b>Applications</b>	ELISA, LF, WB
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements
<b>Size</b>	1 mg
<b>Buffer</b>	10mM PBS, pH 7.4
<b>Preservative</b>	None
<b>Storage</b>	Store at -20°C to -80°C

### BACKGROUND

<b>Keywords</b>	SLA/LP; Soluble Liver Antigen/Liver Pancreas antigen; Autoimmune hepatitis; AIH
-----------------	---