



# Biotinylated SARS-CoV-2 (COVID-19) Spike S2 (BA.1.1) [His, Avi] (DAG-WT817)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Mutations N764K, D796Y, N856K, Q954H, N969K, L981F, F817P, A892P, A899P, A942P, K986P, V987P
<b>Conjugate</b>	His
<b>Applications</b>	Immunoassays
<b>Molecular Weight</b>	62 kDa
<b>Format</b>	Lyophilized
<b>Size</b>	25 µg
<b>Buffer</b>	Lyophilized from PBS, pH 7.4, trehalose
<b>Preservative</b>	None
<b>Storage</b>	Store at -20°C to -80°C

## BACKGROUND

<b>Introduction</b>	SARS-CoV-2 is the virus responsible for COVID-19. It was previously known as 2019 Novel coronavirus (2019-nCoV), is a positive-sense single-stranded RNA virus. SARS-CoV-2 has close genetic similarity to bat coronaviruses, from which it likely originated. An intermediate animal reservoir such as a pangolin is also thought to be involved in its introduction to humans. From a taxonomic perspective, SARS-CoV-2 is classified as a strain of the species severe acute respiratory syndrome-related coronavirus (SARSr-CoV).
<b>Keywords</b>	SARS-CoV-2; Coronavirus; Omicron; SARS-CoV-2 spike protein