



Lactofen [BSA] (DAG-WT837B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Lactofen is a dark brown to tan solid. Insoluble in water. Used as an herbicide.
Keywords	LACTOFEN; 77501-63-4; 2-Ethoxy-1-methyl-2-oxoethyl 5-(2-chloro-4-trifluoromethylphenoxy)-2-nitrobenzoate; (1-ethoxy-1-oxopropan-2-yl) 5-[2-chloro-4-(trifluoromethyl)phenoxy]-2-nitrobenzoate; Cobra