



Repaglinide [BSA] (DAG-WT841B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Repaglinide is a benzoic acid derivative that stimulates insulin secretion from the pancreas and is used in the therapy of type 2 diabetes. Repaglinide has been linked to rare instances of clinically apparent acute liver injury.
Keywords	Repaglinide; 135062-02-1; Prandin; NovoNorm; GlucoNorm