



## Propranolol [OVA] (DAG-WT842O)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Conjugate	OVA
Applications	ELISA, LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C

### BACKGROUND

Introduction	Propranolol is a nonselective beta-adrenergic receptor blocker (beta-blocker) that is widely used for the therapy of hypertension, cardiac arrhythmias, angina pectoris and hyperthyroidism. Propranolol has yet to be convincingly associated with clinically apparent liver injury and is often used in patients with liver disease and cirrhosis.
Keywords	1 Isopropylamino 3 1 Naphthyloxy3Propanol Hydrochloride; Inderal; Inderol; Propranolol HCl; Propranolol Hydrochloride