



## Spironolactone [OVA] (DAG-WT8430)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Conjugate	OVA
Applications	ELISA, LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C

### BACKGROUND

Introduction	Spironolactone is an aldosterone receptor antagonist and potassium-sparing diuretic widely used in the therapy of edema, particularly in patients with cirrhosis in which hyperaldosteronism appears to play a major role. Spironolactone has been linked to rare cases of clinically apparent drug induced liver disease.
Keywords	spironolactone; 52-01-7; Aldactone; Spirolactone; Verospiron