



Disodium 5'-ribonucleotide [OVA] (DAG-WT846O)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	OVA
Applications	ELISA, LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Disodium ribonucleotides is a mixture of disodium inosinate (IMP) and disodium guanylate (GMP) and is often used where a food already contains natural glutamates (as in meat extract) or added monosodium glutamate (MSG). It is primarily used in flavored noodles, snack foods, chips, crackers, sauces and fast foods.
Keywords	Disodium 5'-ribonucleotide; 80702-47-2; tetrasodium;[(2R,3S,4R,5R)-5-(2-amino-6-oxo-1H-purin-9-yl)-3,4-dihydroxyoxolan-2-yl]methyl phosphate;[(2R,3S,4R,5R)-3,4-dihydroxy-5-(6-oxo-1H-purin-9-yl)oxolan-2-yl]methyl phosphate; Disodium5'-ribonucleotide; Disodium 5 Ribonucleotide