



Lysergic acid [OVA] (DAG-WT8490)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	OVA
Applications	ELISA, LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Lysergic acid is an ergoline alkaloid comprising 6-methylergoline having additional unsaturation at the 9,10-position and a carboxy group at the 8-position. It derives from a hydride of a 6-methylergoline.
Keywords	LYSERGIC ACID; 82-58-6; Ergoline Acid; UNII-ITO20DAO7J; ITO20DAO7J