



## Imatinib [OVA] (DAG-WT8500)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Conjugate	OVA
Applications	ELISA, LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C

### BACKGROUND

Introduction	Imatinib is specific tyrosine kinase receptor inhibitor that is used in the therapy of Philadelphia chromosome-positive chronic myelogenous leukemia and gastrointestinal stromal tumors, both of which are marked by an abnormal, constitutively expressed tyrosine kinase that causes unregulated cell growth. Imatinib therapy is associated with transient elevations in serum aminotransferase levels and rare instances of clinically apparent acute liver injury which can be severe and sometimes fatal.
Keywords	Imatinib; 152459-95-5; STI571; Glivec; Imatinib (STI571)