



Tiagabine [BSA] (DAG-WZ011A)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Tiagabine is an anti-convulsive medication. It is also used in the treatment for panic disorder as are a few other anticonvulsants. Though the exact mechanism by which tiagabine exerts its effect on the human body is unknown, it does appear to operate as a selective GABA reuptake inhibitor.
Specificity	N/A
Species	N/A
Conjugate	BSA
Applications	ELISA, LFA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg, 5 mg
Buffer	0.01M PBS, pH 7.2
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	(-)-(R)-1-(4,4-Bis(3-methyl-2-thienyl)-3-butenyl)nipecotic acid; (3R)-1-[4,4-bis(3-methyl-2-thienyl)but-3-en-1-yl]piperidine-3-carboxylic acid; Gabitril; Tiagabinum; Tiagabina; (R)-Tiagabine
-----------------	--