



## Flecainide [BSA] (DAG-WZ016A)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|                            |  |
|----------------------------|--|
| <b>Antigen Description</b> | Flecainide is a potent anti-arrhythmia agent, effective in a wide range of ventricular and atrial arrhythmias and tachycardias. Paradoxically, however, in myocardial infarct patients with either symptomatic or asymptomatic arrhythmia, flecainide exacerbates the arrhythmia and is not recommended for use in these patients. |
| <b>Specificity</b>         | N/A  |
| <b>Species</b>             | N/A  |
| <b>Conjugate</b>           | BSA  |
| <b>Applications</b>        | ELISA, LFA   |
| <b>Format</b>              | Liquid   |
| <b>Concentration</b>       | Batch dependent - please inquire should you have specific requirements.  |
| <b>Size</b>                | 1 mg, 5 mg   |
| <b>Buffer</b>              | 0.01M PBS, pH 7.2  |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.  |

### BACKGROUND

**Keywords** (±)-flecainide; Flecaïne; N-(2-Piperidinylmethyl)-2,5-bis(2,2,2-trifluoroethoxy)benzamide