



Ciclosporin [BSA] (DAG-WZ017A)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Ciclosporin is an immunosuppressant medication and natural product. It is taken by mouth or by injection into a vein for rheumatoid arthritis, psoriasis, Crohn's disease, nephrotic syndrome, and in organ transplants to prevent rejection.
Specificity	N/A
Species	N/A
Conjugate	BSA
Applications	ELISA, LFA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH 7.2
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	CsA; CyA; Cyclosporin; Cyclosporin A; Cyclosporine
-----------------	--