



## Sirolimus [KLH] (DAG-WZ018K)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	10mM PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Wet ice

### BACKGROUND

Introduction	Sirolimus is used in combination with other medications to prevent rejection of kidney transplants. Sirolimus is in a class of medications called immunosuppressants. It works by suppressing the body's immune system.
--------------	---