



SN38 [OVA] (DAG-WZ1009O)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	OVA
Applications	ELISA, LF
Format	Liquid
Size	1 mg
Buffer	0.01M PBS, pH 7.2
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	SN-38 is an antineoplastic drug. It is the active metabolite of irinotecan but has 1000 times more activity than irinotecan itself. In vitro cytotoxicity assays show that the potency of SN-38 relative to irinotecan varies from 2- to 2000-fold.
Keywords	7-Ethyl-10-hydroxycamptothecin; SN 38 lactone