



Berberine[BSA] (DAG-WZ1018)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Berberine is an alkaloid from <i>Hydrastis canadensis</i> L. , Berberidaceae. It is also found in many other plants. It is relatively toxic parenterally, but has been used orally for various parasitic and fungal infections and as antidiarrheal.
Specificity	N/A
Species	N/A
Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg, 5 mg
Buffer	0.01M PBS, pH 7.2
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Berberine; 2086-83-1; Umbellatine; Berberin; UNII-0I8Y3P32UF; 0I8Y3P32UF
-----------------	--