



Capsaicin[BSA] (DAG-WZ1020)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Capsaicin is most often used as a topical analgesic and exists in many formulations of cream, liquid, and patch preparations of various strengths; however, it may also be found in some dietary supplements. Capsaicin is a naturally-occurring botanical irritant in chili peppers, synthetically derived for pharmaceutical formulations.
Specificity	N/A
Species	N/A
Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg, 5 mg
Buffer	0.01M PBS, pH 7.2
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	(E)-8-Methyl-N-vanillyl-6-nonenamide; Isodecenoic acid vanillylamide; trans-8-Methyl-N-vanillyl-6-nonenamide
-----------------	--