



Clindamycin[BSA] (DAG-WZ1024)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Clindamycin is a semisynthetic lincosamide antibiotic that has largely replaced lincomycin due to an improved side effect profile. Clindamycin inhibits bacterial protein synthesis by binding to bacterial 50S ribosomal subunits. It may be bacteriostatic or bactericidal depending on the organism and drug concentration.
Specificity	N/A
Species	N/A
Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg, 5 mg
Buffer	0.01M PBS, pH 7.2
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Clindamycin; Clinacin
-----------------	-----------------------