



## Cyproheptadine[BSA] (DAG-WZ1027)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Cyproheptadine is an antihistamine used to relieve allergy symptoms such as watery eyes, runny nose, itching eyes/nose, sneezing, hives, and itching. It works by blocking a certain natural substance (histamine) that your body makes during an allergic reaction. This medication also blocks another natural substance in your body (serotonin).
<b>Specificity</b>	N/A
<b>Species</b>	N/A
<b>Conjugate</b>	BSA
<b>Applications</b>	ELISA, LFIA
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg, 5 mg
<b>Buffer</b>	0.01M PBS, pH 7.2
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

**Keywords** Axoprol; Ciproheptadina; Cyproheptadin; Cyproheptadine; Cyproheptadinum; Dihexazin