



Diaveridine[BSA] (DAG-WZ1028)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Diaveridine is an aminopyrimidine in which the pyrimidine ring carries amino substituents at C-2 and C-4 and a 3,4-dimethoxybenzyl group at C-5. A folic acid antagonist, it is used as a synergist with sulfonamides against the parasitic Eimeria species. It has a role as an antiparasitic agent.
Specificity	N/A
Species	N/A
Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg, 5 mg
Buffer	0.01M PBS, pH 7.2
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords 5355-16-8; Diaveridin; Diaveridinum; Diaveridina