



## Diaveridine [OVA] (DAG-WZ10280)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Conjugate	OVA
Applications	ELISA, LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	10mM PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Wet ice

### BACKGROUND

Introduction	Diaveridine is an aminopyrimidine in which the pyrimidine ring carries amino substituents at C-2 and C-4 and a 3,4-dimethoxybenzyl group at C-5. A folic acid antagonist, it is used as a synergist with sulfonamides against the parasitic Eimeria species. It has a role as an antiparasitic agent.
--------------	---